

### 便携式多参数重金属分析仪检测标准

序号	重金属离子	测定浓度 (ug/L)	示值误差(%)	重复性 RSD (%)	线性相关系数 r	检测限 (ug/L)
1	铅	10	+0.1	2.6	0.9985	0.1
		50	-2.8	4.6		
2	锌	10	+0.6	4.0	0.9991	0.4
		1000	-1.2	4.8		
3	铜	10	+0.9	4.1	0.9948	0.3
		1000	+2.8	4.5		
4	镉	10	+6.0	3.6	0.9984	0.3
		100	-2.3	3.8		
5	铬 (VI)	25	+2.0	2.8	0.9986	0.9
6	砷 (III)	50	+0.3	4.3	0.9997	0.4
7	汞	10	-0.2	3.4	0.9987	0.1
8	镍	15	0.0	3.7	0.9988	0.2
9	锰	25	+0.8	3.5	0.9985	1.0
10	铊	5	-0.6	4.0	0.9992	0.6
11	钴	10	+0.6	3.7	0.9987	0.2
12	铁	15	+0.7	3.2	0.9992	0.7
13	锑	20	+1.4	2.8	0.9983	1.9